

兔抗 NUP58 多克隆抗体

中文名称：兔抗 NUP58 多克隆抗体

英文名称：Anti-NUP58 rabbit polyclonal antibody

别名：nucleoporin 58; NUP45; NUPL1; PRO2463

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NUP58

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the nucleoporin family that shares 87% sequence identity with rat nucleoporin p58. The protein is localized to the nuclear rim and is a component of the nuclear pore complex (NPC). All molecules entering or leaving the nucleus either diffuse through or are actively transported by the NPC. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Applications:	ELISA, IHC
Name of antibody:	NUP58
Immunogen:	Fusion protein of human NUP58
Full name:	nucleoporin 58
Synonyms:	NUP45; NUPL1; PRO2463

SwissProt:	Q9BVL2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-300



