

OIT3 抗原（重组蛋白）

中文名称： OIT3 抗原（重组蛋白）

英文名称： OIT3 Antigen (Recombinant Protein)

别名： oncoprotein induced transcript 3; LZP

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 122-321 amino acids of human OIT3

技术规格

Full name:	oncoprotein induced transcript 3
Synonyms:	LZP
Swissprot:	Q8WWZ8
Gene Accession:	BC126440
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in hepatocellular function and development.

mlbio 梅联生物
Good elisakit producers