

OLA1 抗原(重组蛋白)

中文名称: OLA1 抗原(重组蛋白)

英文名称: OLA1 Antigen (Recombinant Protein)

别名: Obg like ATPase 1; DOC45; GBP45; GTBP9; GTPBP9; PTD004

储 存: 冷冻(-20℃)

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 197-396 amino acids of human OLA1

技术规格

Full name:	Obg like ATPase 1
Synonyms:	DOC45; GBP45; GTBP9; GTPBP9; PTD004
Swissprot:	Q9NTK5
Gene Accession:	BC012842
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the GTPase protein family. The encoded protein interacts with breast cancer-associated gene 1 (BRCA1) and BRCA1-associated RING domain protein (BARD1), and is involved in centrosome regulation. Overexpression of this gene has been observed in multiple types of cancer and may be a ssociated with poor survival. Pseudogenes of this gene have been defined on chromosomes 17 and 22.