

兔抗 RPS6KB1 (Phospho-Ser424)多克隆抗体

中文名称：兔抗 RPS6KB1 (Phospho-Ser424)多克隆抗体

英文名称：Anti-RPS6KB1 (Phospho-Ser424) rabbit polyclonal antibody

别名：S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：RPS6KB1 (Phospho-Ser424)

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Phosphorylates specifically ribosomal protein S6 in response to insulin or several classes of mitogens. Promotes protein synthesis by phosphorylating PDCD4 at 'Ser-67' and targeting it for degradation.
Applications:	WB, IHC, IF
Name of antibody:	RPS6KB1 (Phospho-Ser424)
Immunogen:	Synthetic peptide of human RPS6KB1 (Phospho-Ser424)
Full name:	ribosomal protein S6 kinase, 70kDa, polypeptide 1 (Phospho-Ser424)
Synonyms :	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-al

	pha; S6K-beta-1; p70(S6K)-alpha
SwissProt:	P23443
IHC positive control:	Human breast carcinoma
IHC Recommend dilution:	50-100
WB Predicted band size:	70 kDa; 85 kDa
WB Positive control:	MCF cells treated with EGF
WB Recommended dilution:	500-1000
IF Positive control:	Hela cells
IF Recommended dilution	100-200



