

兔抗 NUTF2 多克隆抗体

中文名称：兔抗 NUTF2 多克隆抗体

英文名称：Anti-NUTF2 rabbit polyclonal antibody

别名：nuclear transport factor 2; NTF2; PP15; NTF-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NUTF2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62.
Applications:	ELISA, IHC
Name of antibody:	NUTF2
Immunogen:	Fusion protein of human NUTF2
Full name:	nuclear transport factor 2
Synonyms:	NTF2; PP15; NTF-2
SwissProt:	P61970
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	100-300

