

兔抗 INPP1 多克隆抗体

中文名称：兔抗 INPP1 多克隆抗体

英文名称：Anti-INPP1 rabbit polyclonal antibody

别名：inositol polyphosphate-1-phosphatase

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：INPP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the enzyme inositol polyphosphate-1-phosphatase, one of the enzymes involved in phosphatidylinositol signaling pathways. This enzyme removes the phosphate group at position 1 of the inositol ring from the polyphosphates inositol 1,4-bisphosphate and inositol 1,3,4-trisphosphate.
Applications:	ELISA, WB, IHC
Name of antibody:	INPP1
Immunogen:	Fusion protein of human INPP1
Full name:	inositol polyphosphate-1-phosphatase
SwissProt:	P49441

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	44 kDa
WB Positive control:	231 cell and Human fetal liver tissue lysates
WB Recommended dilution:	500-2000



