

兔抗 INPP4A 多克隆抗体

中文名称：兔抗 INPP4A 多克隆抗体

英文名称：Anti-INPP4A rabbit polyclonal antibody

别名：inositol polyphosphate-4-phosphatase type I A; INPP4; TVAS1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：INPP4A

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an Mg ⁺⁺ independent enzyme that hydrolyzes the 4-position phosphate from the inositol ring of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate, and inositol 3,4-bisphosphate. Multiple transcript variants encoding distinct isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	INPP4A
Immunogen:	Fusion protein of human INPP4A
Full name:	inositol polyphosphate-4-phosphatase type I A
Synonyms:	INPP4; TVAS1
SwissProt:	Q96PE3

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer and Human tonsil
IHC Recommend dilution:	100-300



