

兔抗 RPS6KB1(Phospho-Thr389/412) 多克隆抗体

中文名称：兔抗 RPS6KB1(Phospho-Thr389/412) 多克隆抗体

英文名称： Anti-RPS6KB1(Phospho-Thr389/412) rabbit polyclonal antibody

别名： S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha

相关类别： 一抗

储存： 冷冻（-20℃）避光

宿主： Rabbit

抗原： RPS6KB1(Phospho-Thr389/412)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene encodes a member of the ribosomal S6 kinase family of serine/threonine kinases. The encoded protein responds to mTOR (mammalian target of rapamycin) signaling to promote protein synthesis, cell growth

	<p>wth, and cell proliferation. Activity of this gene has been associated with human cancer. Alternatively spliced transcript variants have been observed. The use of alternative translation start sites results in isoforms with longer or shorter N-termini which may differ in their subcellular localizations.</p>
Applications:	WB
Name of antibody:	RPS6KB1(Phospho-Thr389/412)
Immunogen:	Peptide sequence around phosphorylation site of threonine389/412(G-F-T(p)-Y-V) derived from Human p70 S6 Kinase .
Full name:	ribosomal protein S6 kinase, 70kDa, polypeptide 1
Synonyms :	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6K; p70-alpha; S6K-beta-1; p70(S6K)-alpha
SwissProt:	P23443
WB Predicted band size:	59 kDa
WB Positive control:	Jurkat cells lysates treated with Insulin

WB Recommended dilution: 500-1000

