

兔抗 NXPH3 多克隆抗体

中文名称: 兔抗 NXPH3 多克隆抗体

英文名称: Anti-NXPH3 rabbit polyclonal antibody

别名: NPH3

储存: 冷冻 (-20℃) 避光

抗原: NXPH3

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	May be signaling molecules that resemble neuropeptides. Ligand for alpha-neurexins By similarity.
Applications:	WB
Name of antibody:	NXPH3
Immunogen:	Synthesized peptide derived from internal of human NXPH3.
Full name:	neurexophilin 3
Synonyms :	NPH3
SwissProt:	O95157
WB Predicted band size:	28 kDa
WB Positive control:	HT-29 cells lysate
WB Recommended dilution:	500-3000

