

## 兔抗 USP21 多克隆抗体

中文名称：兔抗 USP21 多克隆抗体

英文名称： Anti-USP21 rabbit polyclonal antibody

别名： USP16; USP23

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： USP21

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the C19 peptidase family, also known as family 2 of ubiquitin carboxy-terminal hydrolases. The encoded protein cleaves ubiquitin from ubiquitinated proteins for recycling in intracellular protein degradation. The encoded protein is also able to release NEDD8, a ubiquitin-like protein, from NEDD8-conjugated proteins. This gene has been referred to as USP16 and USP23 but is now known as USP21. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	USP21

<b>Immunogen:</b>	Fusion protein of human USP21
<b>Full name:</b>	ubiquitin specific peptidase 21
<b>Synonyms:</b>	USP16; USP23
<b>SwissProt:</b>	Q9UK80
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	40-200

