

兔抗 NXP4 多克隆抗体

中文名称: 兔抗 NXP4 多克隆抗体

英文名称: Anti-NXP4 rabbit polyclonal antibody

别名: NPH4

储存: 冷冻 (-20℃) 避光

抗原: NXP4

宿主: Rabbit

反应种属: Human Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	May be signaling molecules that resemble neuropeptides and that act by binding to alpha-neurexins and possibly other receptors Potential.
Applications:	WB, IHC
Name of antibody:	NXP4
Immunogen:	Synthesized peptide derived from internal of human NXP4 .
Full name:	neurexophilin 4
Synonyms :	NPH4
SwissProt:	O95158
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

WB Predicted band size:	33 kDa
WB Positive control:	COLO cells and Huvec cells lysates
WB Recommended dilution:	500-3000



