

兔抗 INPP5K 多克隆抗体

中文名称：兔抗 INPP5K 多克隆抗体

英文名称：Anti-INPP5K rabbit polyclonal antibody

别名：inositol polyphosphate-5-phosphatase K; PPS; SKIP; MDCCAID

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：INPP5K

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein with 5-phosphatase activity toward polyphosphate inositol. The protein localizes to the cytosol in regions lacking actin stress fibers. It is thought that this protein may negatively regulate the actin cytoskeleton. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	INPP5K
Immunogen:	Fusion protein of human INPP5K
Full name:	inositol polyphosphate-5-phosphatase K
Synonyms:	PPS; SKIP; MDCCAID

SwissProt:	Q9BT40
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human lung cancer
IHC Recommend dilution:	50-300



