

兔抗 USP24 多克隆抗体

中文名称：兔抗 USP24 多克隆抗体

英文名称： Anti-USP24 rabbit polyclonal antibody

别名： ubiquitin specific peptidase 24

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： USP24

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Modification of cellular proteins by ubiquitin is an essential regulatory mechanism controlled by the coordinated action of multiple ubiquitin-conjugating and deubiquitinating enzymes. USP24 belongs to a large family of cysteine proteases that function as deubiquitinating enzymes (Quesada et al., 2004 [PubMed 14715245]).
Applications:	ELISA, IHC
Name of antibody:	USP24
Immunogen:	Synthetic peptide of human USP24
Full name:	ubiquitin specific peptidase 24

SwissProt:	Q9UPU5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100

