

兔抗 RPS11 多克隆抗体

中文名称：兔抗 RPS11 多克隆抗体

英文名称： Anti-RPS11 rabbit polyclonal antibody

别名： ribosomal protein S11; S11

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPS11

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the S17 P family of ribosomal proteins that is a component of the 40S subunit. This gene is co-transcribed with the small nucleolar RNA gene U35B, which is located in the third intron. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed throughout the genome.

	ome.
Applications:	ELISA, WB, IHC
Name of antibody:	RPS11
Immunogen:	Fusion protein of human RPS11
Full name:	ribosomal protein S11
Synonyms:	S11
SwissProt:	P62280
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human prostate cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	18 kDa
WB Positive control:	Jurkat, NIH/3T3, Raji, 231 and A172 cell lysates
WB Recommended dilution:	500-2000



