

兔抗 RPS12 多克隆抗体

中文名称：兔抗 RPS12 多克隆抗体

英文名称： Anti-RPS12 rabbit polyclonal antibody

别名： ribosomal protein S12; S12

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPS12

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S12E family of ribosomal proteins. It is located in the cytoplasm. Increased expression of this gene in colorectal cancers compared to matched normal colonic mucosa has been observed. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene d

	ispersed through the genome.
Applications:	ELISA, WB, IHC
Name of antibody:	RPS12
Immunogen:	Fusion protein of human RPS12
Full name:	ribosomal protein S12
Synonyms:	S12
SwissProt:	P25398
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	15 kDa
WB Positive control:	Mouse Pancreas tissue lysate
WB Recommended dilution:	500-2000



