

兔抗 USP30 多克隆抗体

中文名称：兔抗 USP30 多克隆抗体

英文名称：Anti-USP30 rabbit polyclonal antibody

别名：ubiquitin specific peptidase 30

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP30

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | USP30, a member of the ubiquitin-specific protease family, is a novel mitochondrial deubiquitinating (DUB) enzyme. |
| Applications: | ELISA, WB |
| Name of antibody: | USP30 |
| Immunogen: | Fusion protein of human USP30 |
| Full name: | ubiquitin specific peptidase 30 |
| SwissProt: | Q70CQ3 |
| ELISA Recommended dilution: | 5000-10000 |
| WB Predicted band size: | 59 kDa |
| WB Positive control: | HepG2 and Hela cell lysates |

WB Recommended dilution: 500-2000

