

## 兔抗 INSC 多克隆抗体

- 中文名称： 兔抗 INSC 多克隆抗体
- 英文名称： Anti-INSC rabbit polyclonal antibody
- 别名： INSC, spindle orientation adaptor protein
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 抗原： INSC
- 反应种属： Human, Mouse
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	In <i>Drosophila</i> , neuroblasts divide asymmetrically into another neuroblast at the apical side and a smaller ganglion mother cell on the basal side. Cell polarization is precisely regulated by 2 apically localized multiprotein signaling complexes that are tethered by Inscuteable, which regulates their apical localization (Izaki et al., 2006 [PubMed 16458856]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	INSC
<b>Immunogen:</b>	Synthetic peptide of human INSC
<b>Full name:</b>	INSC, spindle orientation adaptor protein

<b>SwissProt:</b>	Q1MX18
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human liver cancer
<b>IHC Recommend dilution:</b>	50-300



