

## 兔抗 OAS3 多克隆抗体

中文名称：兔抗 OAS3 多克隆抗体

英文名称：Anti-OAS3 rabbit polyclonal antibody

别名：2'-5'-oligoadenylate synthetase 3; p100; p100OAS

相关类别：一抗

抗原：OAS3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an enzyme included in the 2', 5' oligo adenylate synthase family. This enzyme is induced by interferons and catalyzes the 2', 5' oligomers of adenosine in order to bind and activate RNase L. This enzyme family plays a significant role in the inhibition of cellular protein synthesis and viral infection resistance.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	OAS3
<b>Immunogen:</b>	Synthetic peptide of human OAS3
<b>Full name:</b>	2'-5'-oligoadenylate synthetase 3
<b>Synonyms:</b>	p100; p100OAS
<b>SwissProt:</b>	Q9Y6K5
<b>ELISA Recommended dilution:</b>	5000-10000

IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-10

