

兔抗 INSM1 多克隆抗体

中文名称: 兔抗 INSM1 多克隆抗体

英文名称: Anti-INSM1 rabbit polyclonal antibody

别名: INSM transcriptional repressor 1; IA1; IA-1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: INSM1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Insulinoma-associated 1 (INSM1) gene is intronless and encodes a protein containing both a zinc finger DNA-binding domain and a putative prohormone domain. This gene is a sensitive marker for neuroendocrine differentiation of human lung tumors.
Applications:	ELISA, IHC
Name of antibody:	INSM1
Immunogen:	Synthetic peptide of human INSM1
Full name:	INSM transcriptional repressor 1
Synonyms:	IA1; IA-1

SwissProt:	Q01101
ELISA Recommended dilution :	5000-10000
IHC positive control:	Human esophagus cancer and Human breast cancer
IHC Recommend dilution:	50-300

