

兔抗 OASL 多克隆抗体

中文名称：兔抗 OASL 多克隆抗体

英文名称： Anti-OASL rabbit polyclonal antibody

别名： OASLd, TRIP14, TRIP-14, p59OASL, p59 OASL, p59-OASL

抗原： OASL

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	OASL encodes an interferon-induced protein that is a member of the 2-5A synthetase family. 2-prime,5-prime oligoadenylates (2-5As) bind to and activate RNase L, resulting in general RNA degradation and consequent inhibition of protein synthesis. 2-5As are produced by a well-conserved family of interferon-induced enzymes, the 2-5A synthetases or OASs
Applications:	ELISA, IHC
Name of antibody:	OASL
Immunogen:	Fusion protein of human OASL
Full name:	2'-5'-oligoadenylate synthetase-like
Synonyms :	OASLd, TRIP14, TRIP-14, p59OASL, p59 OASL, p59-OASL
SwissProt:	Q15646

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	25-100

