

兔抗 RPS15A 多克隆抗体

中文名称：兔抗 RPS15A 多克隆抗体

英文名称： Anti-RPS15A rabbit polyclonal antibody

别名： S15a

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPS15A

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S8P family of ribosomal proteins. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

Applications:	ELISA, WB, IHC
Name of antibody:	RPS15A
Immunogen:	Fusion protein of human RPS15A
Full name:	ribosomal protein S15a
Synonyms:	S15a
SwissProt:	P62244
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human colorectal cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	15 kDa
WB Positive control:	Human muscle tissue lysate
WB Recommended dilution:	200-1000



