

## 兔抗 USP38 多克隆抗体

中文名称：兔抗 USP38 多克隆抗体

英文名称：Anti-USP38 rabbit polyclonal antibody

别名：HP43.8KD

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：USP38

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	Deubiquitinating enzyme exhibiting a preference towards 'Lys-63'-linked Ubiquitin chains.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	USP38
<b>Immunogen:</b>	Synthesized peptide derived from internal of human USP38.
<b>Full name:</b>	ubiquitin specific peptidase 38
<b>Synonyms :</b>	HP43.8KD
<b>SwissProt:</b>	Q8NB14
<b>WB Predicted band size:</b>	117 kDa
<b>WB Positive control:</b>	293 cells lysate

**WB Recommended dilution:** 500-3000

