

兔抗 OCLN 多克隆抗体

中文名称：兔抗 OCLN 多克隆抗体

英文名称： Anti-OCLN rabbit polyclonal antibody

别名： BLCPMG

抗原： OCLN

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an integral membrane protein that is required for cytokine-induced regulation of the tight junction paracellular permeability barrier. Mutations in this gene are thought to be a cause of band-like calcification with simplified gyration and polymicrogyria (BLC-PMG), an autosomal recessive neurologic disorder that is also known as pseudo-TORCH syndrome. Alternative splicing results in multiple transcript variants. A related pseudogene is present 1.5 Mb downstream on the q arm of chromosome 5.
Applications:	ELISA, IHC
Name of antibody:	OCLN
Immunogen:	Fusion protein of human OCLN

Full name:	Occludin
Synonyms :	BLCPMG
SwissProt:	Q16625
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer and Human lung cancer
IHC Recommend dilution:	15-50

