

## 兔抗 RPS16 多克隆抗体

中文名称：兔抗 RPS16 多克隆抗体

英文名称： Anti-RPS16 rabbit polyclonal antibody

别名： ribosomal protein S16; S16

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RPS16

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S9P family of ribosomal proteins. It is located in the cytoplasm. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RPS16
<b>Immunogen:</b>	Fusion protein of human RPS16
<b>Full name:</b>	ribosomal protein S16
<b>Synonyms:</b>	S16
<b>SwissProt:</b>	P62249
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	100-200



