

## 兔抗 USP42 多克隆抗体

中文名称：兔抗 USP42 多克隆抗体

英文名称：Anti-USP42 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）避光

宿 主：Rabbit

抗 原：USP42

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Deubiquitinating enzyme which may play an important role during spermatogenesis By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	USP42
<b>Immunogen:</b>	Synthesized peptide derived from internal of human USP42 .
<b>Full name:</b>	ubiquitin specific peptidase 42
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q9H9J4
<b>WB Predicted band size:</b>	145 kDa
<b>WB Positive control:</b>	HT-29 cells lysate
<b>WB Recommended dilution:</b>	500-3000

