

兔抗 INSR(Phospho-Tyr1355) 多克隆抗体

中文名称：兔抗 INSR(Phospho-Tyr1355) 多克隆抗体

英文名称：Anti-INSR(Phospho-Tyr1355) rabbit polyclonal antibody

别名：HHF5; CD220

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：INSR(Phospho-Tyr1355)

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	After removal of the precursor signal peptide, the insulin receptor precursor is post-translationally cleaved into two chains (alpha and beta) that are covalently linked. Binding of insulin to the insulin receptor (INSR) stimulates glucose uptake. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	INSR(Phospho-Tyr1355)
Immunogen:	Peptide sequence around phosphorylation site of tyrosine 1355 (R-S-Y(p)-E-E) derived from Human IR.

Full name:	insulin receptor
Synonyms :	HHF5; CD220
SwissProt:	P06213
WB Predicted band size:	113 kDa
WB Positive control:	293 cells lysates treated with Heatshock
WB Recommended dilution:	500-1000

