

兔抗 RPS20 多克隆抗体

中文名称：兔抗 RPS20 多克隆抗体

英文名称： Anti-RPS20 rabbit polyclonal antibody

别名： S20

抗原： RPS20

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S10P family of ribosomal proteins. It is located in the cytoplasm. This gene is co-transcribed with the small nucleolar RNA gene U54, which is located in its second intron. As is typical for genes encoding ribosomal proteins, there are

	multiple processed pseudogenes of this gene dispersed through the genome. Two transcript variants encoding different isoforms have been identified for this gene.
Applications:	WB, IF
Name of antibody:	RPS20
Immunogen:	Synthesized peptide derived from internal of human RPS20.
Full name:	ribosomal protein S20
Synonyms :	S20
SwissProt:	P60866
WB Predicted band size:	13 kDa
WB Positive control:	Huvec cells and MCF7 cells lysates
WB Recommended dilution:	500-3000
IF positive control:	A549 cells
IF Recommend dilution:	100-500



