

兔抗 USP47 多克隆抗体

中文名称：兔抗 USP47 多克隆抗体

英文名称：Anti-USP47 rabbit polyclonal antibody

别名：TRFP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP47

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Ubiquitin carboxyl-terminal hydrolase 47 is an enzyme that in humans is encoded by the USP47 gene. Ubiquitin-specific protease that specifically deubiquitinates monoubiquitinated DNA polymerase beta (POLB), stabilizing POLB thereby playing a role in base-excision repair (BER). Acts as a regulator of cell growth and genome integrity. May also indirectly regulates CDC25A expression at a transcriptional level.
Applications:	ELISA, IHC
Name of antibody:	USP47
Immunogen:	Fusion protein of human USP47

Full name:	ubiquitin specific peptidase 47
Synonyms :	TRFP
SwissProt:	Q96K76
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	25-100

