

## 兔抗 OGFR 多克隆抗体

中文名称: 兔抗 OGFR 多克隆抗体

英文名称: Anti-OGFR rabbit polyclonal antibody

抗 原: OGFR

储 存: 冷冻 (-20℃) 避光

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

<b>Background:</b>	The protein encoded by this gene is a receptor for opioid growth factor (OGF), also known as [Met(5)]-enkephalin. OGF is a negative regulator of cell proliferation and tissue organization in a variety of processes. The encoded unbound receptor for OGF has been localized to the outer nuclear envelope, where it binds OGF and is translocated into the nucleus. The coding sequence of this gene contains a polymorphic region of 60 nt tandem imperfect repeat units. Several transcripts containing between zero and eight repeat units have been reported.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	OGFR
<b>Immunogen:</b>	Synthesized peptide derived from internal of human OGFR.
<b>Full name:</b>	opioid growth factor receptor
<b>Synonyms :</b>	0

<b>SwissProt:</b>	Q9NZT2
<b>WB Predicted band size:</b>	73 kDa
<b>WB Positive control:</b>	Hela cells and HepG2 cells lysates
<b>WB Recommended dilution:</b>	500-3000

