

兔抗 OIT3 多克隆抗体

中文名称：兔抗 OIT3 多克隆抗体

英文名称：Anti-OIT3 rabbit polyclonal antibody

别名：LZP

抗原：OIT3

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in hepatocellular function and development.
Applications:	ELISA, IHC
Name of antibody:	OIT3
Immunogen:	Fusion protein of human OIT3
Full name:	oncogene induced transcript 3
Synonyms :	LZP
SwissProt:	Q8WWZ8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

