

兔抗 USP53 多克隆抗体

中文名称：兔抗 USP53 多克隆抗体

英文名称：Anti-USP53 rabbit polyclonal antibody

别名：ubiquitin specific peptidase 53

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：USP53

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Tight junction-associated protein that is involved in the survival of auditory hair cells and hearing. Maybe by modulating the barrier properties and mechanical stability of tight junctions (By similarity). Has no peptidase activity (PubMed:14715245).
Applications:	ELISA, IHC
Name of antibody:	USP53
Immunogen:	Fusion protein of human USP53
Full name:	ubiquitin specific peptidase 53
SwissProt:	Q70EK8
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human cervical cancer
IHC Recommend dilution:	100-300

