

兔抗 RPS27L 多克隆抗体

中文名称: 兔抗 RPS27L 多克隆抗体

英文名称: Anti-RPS27L rabbit polyclonal antibody

别名: ribosomal protein S27 like

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: RPS27L

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein sharing 96% amino acid similarity with ribosomal protein S27, which suggests the encoded protein may be a component of the 40S ribosomal subunit.
Applications:	ELISA, WB, IHC
Name of antibody:	RPS27L
Immunogen:	Fusion protein of human RPS27L
Full name:	ribosomal protein S27 like
SwissProt:	Q71UM5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human thyroid ca

	ncer
IHC Recommend dilution:	40-200
WB Predicted band size:	9 kDa
WB Positive control:	HT-29 cell, LOVO cell, 231 cell, Rat lung tissue, A172 cell, HUVEC cell and Mouse adrenal gland tissue lysates
WB Recommended dilution:	200-1000





