

PDRG1 抗原（重组蛋白）

中文名称：PDRG1 抗原（重组蛋白）

英文名称：PDRG1 Antigen (Recombinant Protein)

别名：p53 and DNA damage regulated 1; PDRG; C20orf126

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-133 amino acids of human PDRG1

技术规格

Full name:	p53 and DNA damage regulated 1
Synonyms:	PDRG; C20orf126
Swissprot:	Q9NUG6
Gene Accession:	BC009334
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May play a role in chaperone-mediated protein folding.