

## 兔抗 UTP4 多克隆抗体

中文名称：兔抗 UTP4 多克隆抗体

英文名称：Anti-UTP4 rabbit polyclonal antibody

别名：UTP4, small subunit processome component; NAIC; CIRH1A; CIRHIN; TEX292

相关类别：一抗

抗原：UTP4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a WD40-repeat-containing protein that is localized to the nucleolus. Mutation of this gene causes North American Indian childhood cirrhosis, a severe intrahepatic cholestasis that results in transient neonatal jaundice, and progresses to periportal fibrosis and cirrhosis in childhood and adolescence. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UTP4
<b>Immunogen:</b>	Synthetic peptide of human UTP4
<b>Full name:</b>	UTP4, small subunit processome component

<b>Synonyms:</b>	NAIC; CIRH1A; CIRHIN; TEX292
<b>SwissProt:</b>	Q969X6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	30-150

