

## 兔抗 OLFM4 多克隆抗体

中文名称：兔抗 OLFM4 多克隆抗体

英文名称： Anti-OLFM4 rabbit ployclonal antibody

别 名： GC1; OLM4; OlfD; GW112; hGC-1; hOLfD; UNQ362; bA209J19.1

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： OLFM4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene was originally cloned from human myeloblasts and found to be selectively expressed in inflamed colonic epithelium. This gene encodes a member of the olfactomedin family. The encoded protein is an antiapoptotic factor that promotes tumor growth and is an extracellular matrix glycoprotein that facilitates cell adhesion.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	OLFM4
<b>Immunogen:</b>	Fusion protein of human OLFM4
<b>Full name:</b>	olfactomedin 4
<b>Synonyms :</b>	GC1; OLM4; OlfD; GW112; hGC-1; hOLfD; UNQ362; bA209J1

	9.1
<b>SwissProt:</b>	Q6UX06
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	57 kDa
<b>WB Positive control:</b>	Human fetal liver tissue
<b>WB Recommended dilution:</b>	200-1000

