

兔抗 UTP18 多克隆抗体

中文名称：兔抗 UTP18 多克隆抗体

英文名称：Anti-UTP18 rabbit polyclonal antibody

别名：WDR50; CGI-48

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UTP18

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Involved in nucleolar processing of pre-18S ribosomal RNA.
Applications:	ELISA, WB, IHC
Name of antibody:	UTP18
Immunogen:	Synthetic peptide of human UTP18
Full name:	UTP18 small subunit processome component
Synonyms:	WDR50; CGI-48
SwissProt:	Q9Y5J1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer

IHC Recommend dilution:	20-100
WB Predicted band size:	62 kDa
WB Positive control:	293T , Hela , Raji and A172 cell lysates
WB Recommended dilution:	200-1000



