

兔抗 OLFML3 多克隆抗体

中文名称：兔抗 OLFML3 多克隆抗体

英文名称： Anti-OLFML3 rabbit polyclonal antibody

别名： olfactomedin like 3; OLF44; HNOEL-iso

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： OLFML3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Secreted scaffold protein that plays an essential role in dorsoventral patterning during early development. Stabilizes axial formation by restricting chordin (CHRD) activity on the dorsal side. Acts by facilitating the association between the tolloid proteases and their substrate chordin (CHRD), leading to enhance chordin (CHRD) degradation (By similarity). May have matrix-related function involved in placental and embryonic development, or play a similar role in other physiological processes.
Applications:	ELISA, IHC
Name of antibody:	OLFML3

Immunogen:	Fusion protein of human OLFML3
Full name:	olfactomedin like 3
Synonyms:	OLF44; HNOEL-iso
SwissProt:	Q9NRN5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

