

## PDGFB 抗原（多肽片段）

中文名称： PGR 抗原（多肽片段）

英文名称： PGR Antigen

相关类别： 抗原

储 存： 冷冻（-20℃）

### 概述

Synthetic peptide corresponding to a region derived from 433-447 amino acids of human progesterone receptor

### 技术规格

<b>Full name:</b>	Progesterone receptor
<b>Synonyms:</b>	PR, NR3C3
<b>Gene Accession:</b>	NP_001258090
<b>SwissProt:</b>	P06401
<b>Background:</b>	This gene encodes a member of the steroid receptor superfamily. The encoded protein mediates the physiological effects of progesterone, which plays a central role in reproductive events associated with the establishment and maintenance of pregnancy. This gene uses two distinct promoters and translation start sites in the first exon to produce two isoforms, A and B. The two isoforms are identical except for the additional 165 amino acids found in the N-terminus of isoform B and mediate their own response genes and physiologic effects with little overlap.

**mlbio** 梅联生物  
Good elisakit producers