

兔抗 OPRD1 多克隆抗体

中文名称：兔抗 OPRD1 多克隆抗体

英文名称： Anti-OPRD1 rabbit polyclonal antibody

别名： OPRD

抗原： OPRD1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The opioid receptor, delta 1, also known as delta opioid receptor or simply delta receptor, abbreviated DOR, is an opioid receptor that has enkephalins as its endogenous ligands. Pharmacological studies and molecular cloning have led to the identification of three different types of opioid receptor, mu-type, delta-type and kappa-type, also designated MOR-1, DOR-1 and KOR-1, respectively. MOR-1 is a receptor for beta-endorphin, DOR-1 is a receptor for enkephalins, and KOR-1 is a receptor for dynorphins.
Applications:	ELISA, IHC

Name of antibody:	OPRD1
Immunogen:	Synthetic peptide of human OPRD1
Full name:	opioid receptor, delta 1
Synonyms :	OPRD
SwissProt:	P41143
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and Human breast cancer
IHC Recommend dilution:	25-100

