

兔抗 UTS2 多克隆抗体

中文名称：兔抗 UTS2 多克隆抗体

英文名称：Anti-UTS2 rabbit polyclonal antibody

别名：U11, U-11, UCN2

相关类别：一抗

抗原：UTS2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a mature peptide that is an active cyclic heptapeptide absolutely conserved from lamprey to human. The active peptide acts as a vasoconstrictor and is expressed only in brain tissue. Despite the gene family name similarity, this gene is not homologous to urocortin, a member of the sauvagine/corticotropin-releasing factor/urotensin I family. Most of the proprotein is cleaved to make the mature peptide. Transcript variants encoding different preproprotein isoforms have been described for this gene.

Applications:	ELISA, IHC
Name of antibody:	UTS2
Immunogen:	Fusion protein of human UTS2
Full name:	Urotensin 2
Synonyms :	UII, U-II, UCN2
SwissProt:	O95399
ELISA Recommended dilution:	500-5000
IHC positive control:	Human lung cancer and human cervical cancer
IHC Recommend dilution:	50-200

