

兔抗 IPO13 多克隆抗体

中文名称：兔抗 IPO13 多克隆抗体

英文名称：Anti-IPO13 rabbit polyclonal antibody

别名：LGL2; IMP13; KAP13; RANBP13

相关类别：一抗

储存：冷冻（-20℃）

抗原：IPO13

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the importin-beta family of nuclear transport proteins. The encoded protein mediates the import of specific cargo proteins from the cytoplasm to the nucleus and is dependent on the Ras-related nuclear protein-GTPase system. The encoded protein is also involved in nuclear export of the eukaryotic translation initiation factor 1A.
Applications:	ELISA, IHC
Name of antibody:	IPO13

Immunogen:	Fusion protein of human IPO13
Full name:	importin 13
Synonyms :	LGL2; IMP13; KAP13; RANBP13
SwissProt:	O94829
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human colon cancer
IHC Recommend dilution:	50-200

