

兔抗 UTY 多克隆抗体

中文名称: 兔抗 UTY 多克隆抗体

英文名称: Anti-UTY rabbit polyclonal antibody

别名: UTY1; KDM6AL

抗原: UTY

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a protein containing tetratricopeptide repeats which are thought to be involved in protein-protein interactions. The encoded protein is also a minor histocompatibility antigen which may induce graft rejection of male stem cell grafts. A large number of alternatively spliced transcripts have been observed for this gene, but the full length nature of some of these variants has not been determined.

Applications:	ELISA, IHC
Name of antibody:	UTY
Immunogen:	Full length fusion protein
Full name:	ubiquitously transcribed tetratricopeptide repeat c ontaining, Y-linked
Synonyms:	UTY1; KDM6AL
SwissProt:	O14607
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and human thyroid cancer
IHC Recommend dilution:	25-100



