


## 兔抗 UXS1 多克隆抗体

- 中文名称：兔抗 UXS1 多克隆抗体
- 英文名称：Anti-UXS1 rabbit polyclonal antibody
- 别名：UDP-glucuronate decarboxylase 1; UGD; SDR6E1
- 储存：冷冻（-20℃）
- 抗原：UXS1
- 宿主：Rabbit
- 反应种属：Human, Mouse, Rat
- 相关类别：一抗
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an enzyme found in the perinuclear Golgi which catalyzes the synthesis of UDP-xylose used in glycosaminoglycan (GAG) synthesis on proteoglycans. The GAG chains are covalently attached to proteoglycans which participate in signaling pathways during development. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	UXS1
<b>Immunogen:</b>	Fusion protein of human UXS1
<b>Full name:</b>	UDP-glucuronate decarboxylase 1
<b>Synonyms:</b>	UGD; SDR6E1
<b>SwissProt:</b>	Q8NBZ7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human tonsil
<b>IHC Recommend dilution:</b>	50-300



**mlbio** 梅联生物  
Good elisakit producers

