

兔抗 IQCC 多克隆抗体

中文名称: 兔抗 IQCC 多克隆抗体

英文名称: Anti-IQCC rabbit polyclonal antibody

别名: IQ motif containing C

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: IQCC

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, WB, IHC
Name of antibody:	IQCC
Immunogen:	Fusion protein of human IQCC
Full name:	IQ motif containing C
SwissProt:	Q4KMZ1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	150-300
WB Predicted band size:	53 kDa

WB Positive control:	LOVO cell lysate
WB Recommended dilution:	1000-5000



