

兔抗 UXT 多克隆抗体

中文名称：兔抗 UXT 多克隆抗体

英文名称：Anti-UXT rabbit polyclonal antibody

别名：STAP1; ART-27

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：UXT

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene functions as a cofactor that modulates androgen receptor-dependent transcription, and also plays a critical role in tumor necrosis factor-induced apoptosis. Expression of this gene may play a role in tumorigenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, IHC
Name of antibody:	UXT
Immunogen:	Full length fusion protein
Full name:	ubiquitously-expressed, prefoldin-like chaperone

Synonyms :	STAP1; ART-27
SwissProt:	Q9UBK9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer
IHC Recommend dilution:	20-100

