

## 兔抗 RSAD1 多克隆抗体

中文名称: 兔抗 RSAD1 多克隆抗体

英文名称: Anti-RSAD1 rabbit polyclonal antibody

抗原: RSAD1

储存: 冷冻 (-20°C) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	May be involved in porphyrin cofactor biosynthesis By similarity.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	RSAD1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human RSAD1.
<b>Full name:</b>	radical S-adenosyl methionine domain containing 1
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q9HA92
<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	49 kDa

<b>WB Positive control:</b>	Hela cells lysate
<b>WB Recommended dilution:</b>	500-3000



