

兔抗 IQSEC1 多克隆抗体

中文名称：兔抗 IQSEC1 多克隆抗体

英文名称：Anti-IQSEC1 rabbit polyclonal antibody

别名：IQ motif and Sec7 domain 1; BRAG2; GEP100; ARFGEP100; ARF-GEP100

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IQSEC1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Guanine nucleotide exchange factor for ARF1 and ARF6 (PubMed:24058294). Guanine nucleotide exchange factor activity is enhanced by lipid binding (PubMed:24058294). Accelerates GTP binding by ARFs of all three classes. Guanine nucleotide exchange protein for ARF6, mediating internalisation of beta-1 integrin.
Applications:	ELISA, IHC
Name of antibody:	IQSEC1
Immunogen:	Fusion protein of human IQSEC1
Full name:	IQ motif and Sec7 domain 1

Synonyms:	BRAG2; GEP100; ARFGEP100; ARF-GEP100
SwissProt:	Q6DN90
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human gastric cancer
IHC Recommend dilution:	150-300



