

## 兔抗 RSPH14 多克隆抗体

中文名称：兔抗 RSPH14 多克隆抗体

英文名称： Anti-RSPH14 rabbit polyclonal antibody

别名： radial spoke head 14 homolog; RTDR1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RSPH14

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with no known function but with slight similarity to a yeast vacuolar protein. The gene is located in a region deleted in pediatric rhabdoid tumors of the brain, kidney and soft tissues, but mutations in this gene have not been associated with the disease.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RSPH14
<b>Immunogen:</b>	Fusion protein of human RSPH14
<b>Full name:</b>	radial spoke head 14 homolog
<b>Synonyms:</b>	RTDR1

<b>SwissProt:</b>	Q9UHP6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100

