

## 兔抗 RSU1 多克隆抗体

中文名称：兔抗 RSU1 多克隆抗体

英文名称： Anti-RSU1 rabbit polyclonal antibody

别名： RSP-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RSU1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is involved in the Ras signal transduction pathway, growth inhibition, and nerve-growth factor induced differentiation processes, as determined in mouse and human cell line studies. In mouse, the encoded protein was initially isolated based on its ability to inhibit v-Ras transformation. Multiple alternatively spliced transcript variants for this gene have been reported; one of these variants was found only in glioma tumors.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RSU1

<b>Immunogen:</b>	Fusion protein of human RSU1
<b>Full name:</b>	Ras suppressor protein 1
<b>Synonyms:</b>	RSP-1
<b>SwissProt:</b>	Q15404
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	32 kDa
<b>WB Positive control:</b>	Rat heart tissue, Mouse heart tissue, Rat brain tissue, Mouse brain tissue, HEPG2 cell, Jurkat and Lovo cell lysates
<b>WB Recommended dilution:</b>	500-2000



