

兔抗 VAMP4 多克隆抗体

中文名称: 兔抗 VAMP4 多克隆抗体

英文名称: Anti-VAMP4 rabbit polyclonal antibody

别名: VAMP-4; VAMP24

储存: 冷冻 (-20℃)

抗原: VAMP4

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. This protein may play a role in trans-Golgi network-to-endosome transport.
Applications:	ELISA, WB, IHC
Name of antibody:	VAMP4
Immunogen:	Full length fusion protein

Full name:	vesicle-associated membrane protein 4
Synonyms :	VAMP-4; VAMP24
SwissProt:	O75379
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	16 kDa
WB Positive control:	PC3, TM4, hela and K562 cell
WB Recommended dilution:	500-2000



